



## Carestream CR Mobile Installation Kit

Retrofits mosts vehicles



Pictured here in a Chevy HHR



Can be customized to install in the front or back of any vehicle!

### Mobile CR Installation Kit

This unique kit includes a custom designed fiberglass enclosure that creates a desk top surface. Underneath the enclosure is shock absorbing material and tie down straps that protect the CR Reader from moving. The laptop is mounted and locked on the enclosure for added protection. To power the CR a 1500W Pure Sine Wave Inverter is provided along with a battery isolator, circuit breaker and on/off switch.

A detailed wiring diagram and step by step instructions make installing this kit quick and easy!